





Endoscope Bedside Kit



StageOne™ is an all-in-one Tub, detergent channel cleanser and scope sponge designed for the initial flush and wipe of an endoscope immediately post-procedure. It consists of a unique biodegradable pot, and a sachet containing a universal sponge and detergent concentrate. All you need to do is just add water!

Features:

- Addresses the need of international guidelines
- Ecological tub. Bagasse is an eco-friendly, biodegradable, renewable organic material. In addition it is nontoxic, harmless, healthy and sanitary, and can be recycled and good to the environment. Bagasse is a renewable and sustainable resource and is 100% compostable.
- Easy peel sachet and universal shaped reticulated sponge
- The unique two-level design of the StageOne™ tub allows for either clinical or Olympus guidelines to be followed
- Easy to transport and store
- No water, reducing the weight. StageOne™ is exceptionally light, reducing the burden on transport costs by up-to 30% and the product uses less packing material reducing the waste stream



The StageOne Bedside Cleaning kit addresses all the recent US, World and European Guidelines.

AAMI ST91;2015

To prevent build-up of bioburden, development of biofilms, and drying of secretions, precleaning should take place at the point of use immediately following the procedure

SGNA:

Pre-cleaning should be performed at the point of use, before bioburden has an opportunity to dry and before complete decontamination.

WGO/OMED

Pre-cleaning: Clear gross debris by sucking detergent through the working channel (250 ml).

ESGE


Before the endoscope is detached from the light source and video processor, detergent solution is sucked through the working channel. Wipe down the insertion tube with a soft, disposable cloth/ sponge.

Olympus

Immerse the distal end of the insertion section in the solution. Depress the suction valve on the endoscope and aspirate detergent through the endoscope for 30 seconds.

This equates to about 500ml of detergent solution.

Instructions for Use



Product	Part Number	Description
StageOne™	INOV8-009	Stage-One™ : Endoscope Ultra Pre-Cleaner for the initial flush and wipe of an endoscope. Unique 250ml and 500ml two level pot. Includes integral scope sponge. Pack of 100