

MEDIVATORS®

# ENDO STRATUS®

Irrigation Pump and CO2 Insufflator




MEDIVATORS®

# ENDO STRATUS®

Irrigation Pump and CO<sub>2</sub> Insufflator

Patient Pain and Discomfort CO<sub>2</sub> vs Room Air



CO<sub>2</sub> Insufflation is proven to result in faster procedures, quicker recovery times and more comfortable patients.

## Quicker Patient Turnaround\*

- Improved cecal intubation rates<sup>1</sup>
- Greater small bowel intubation depths<sup>2</sup>
- 38% reduction in nursing attention<sup>4</sup>

## Higher Patient Satisfaction

- Less discomfort, pain and bloating<sup>3,4,1</sup>
- Quicker recovery times<sup>4</sup>
- Decreased gas distension

1. Singh, R., Neo, E., & Nordeen, N. (2012, July). Carbon dioxide insufflation during colonoscopy in deeply sedated patients. *World Journal of Gastroenterology*, 18(25), 3250-3253.

2. Bretthauer, MD, PhD, M. (2010). Turning science into clinical practice – the case of carbon dioxide insufflation. *Endoscopy* 2010, 42(9), 1104-1105. doi:http://dx.doi.org/ 10.1055/s-0030-1255973

3. Bretthauer, MD, PhD, M. (2007, September). Carbon dioxide insufflation improves intubation depth in double-balloon enteroscopy: a randomized, controlled, double-blind trial. *Endoscopy* 2007, 39(9), 1064-1067. doi:DOI 10.1055/s-2007-966990

4. Lynch, MBA, BSN, RN, Hayes BSN, RN, CGRN, Buffum DNSc, APRN, CS. (2012). Insufflation Using CO<sub>2</sub> vs Room Air During Colonoscopy: Comparison of Patient Comfort, Recovery Time, and Nursing Resources [PowerPoint slides]. Veterans Affairs Medical Center San Francisco.

# Simple Change, Significant Improvement

## ENDO STRATUS® CO<sub>2</sub> Insufflator

- Absorbed 150 times faster than room air and is promptly eliminated via the lungs
- Heated CO<sub>2</sub> – 98.6 degrees
- Variable flow rates – low, medium, & high

### Compatible with

- MEDIVATORS® ENDO SMARTCAP® lens cleaning and insufflation tubing
- Tank (D or E) or DISS Wall CO<sub>2</sub>
- All major GI Endoscopes



# Improved Visualization and Patient Comfort


## ENDO STRATUS® Irrigation Pump

- Warm water technology minimizes chances of spastic colon
- Automatic prime button provides instant irrigation upon foot pedal depression
- Irrigate through auxillary water jet or biopsy channel

### Compatible with

- MEDIVATORS® ENDOGATOR® tubing and connectors
- All major GI Endoscopes






## CO<sub>2</sub> Carts

- Conveniently transport and store two CO<sub>2</sub> tanks.
- ENDO STRATUS® Irrigation Pump and CO<sub>2</sub> Insufflation units stack neatly on top platform.
- Non-porous straps secure ENDO STRATUS Irrigation Pump and CO<sub>2</sub> Insufflation units for easy transportation between rooms.
- Basket provides additional storage.

### ENDO STRATUS™ Carts Ordering Information

Part Number	Description
EGA-7030	ENDO STRATUS IV Style Cart
EGA-7032	ENDO STRATUS IV Style Cart with Dual Yoke
EGA-7031	ENDO STRATUS Box Style Cart
EGA-7033	ENDO STRATUS Box Style Cart with Dual Yoke



### ENDO STRATUS® Pumps Ordering Information

Part Number	Description
EGA-500	ENDO STRATUS Irrigation Pump
EGA-501	ENDO STRATUS Insufflator Unit
EGA-SYST	ENDO STRATUS Irrigation Pump & CO <sub>2</sub> Insufflator
EGA-7012	Tank Yoke
EGA-7029	Dual Tank Yoke
EGA-7020	3 Foot Wall Hose (Low Pressure)
EGA-7023	9 Foot Wall Hose (Low Pressure)
EGA-7011	3 Foot Tank Hose (High Pressure)
EGA-7024	9 Foot Tank Hose (High Pressure)

MEDIVATORS®, ENDO STRATUS®, ENDO SMARTCAP® and ENDOGATOR® are registered trademarks of Medivators Inc.



MEDIVATORS Inc. Headquarters  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.800.444.4729  
Fax: +1.763.553.3387

MEDIVATORS BV  
Sorethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte. Ltd.  
1 A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office  
Room 708, 7th Floor Kaiheng Center, Block B  
No. 2 Chaoyangmennei Street  
Dongcheng District, Beijing China 100010  
Tel: +8610.6567.8446  
Fax: +8610.6567.8445

www.medivators.com

© 2014 MEDIVATORS Inc. P/N 50098-583/A