

Eclipse Piston

The Rebirth of Nitinol



Grace Medical
ENT. Redefined.

Shape Memory Zone



Memory Only
Where Needed



Complete Gentle
Closure Around Incus

No Movement
in Vestibule



*De la Cruz Style

Cat.#	Shaft Dia.	Length
469-400	0.6mm	4.00mm
469-425	0.6mm	4.25mm
469-450	0.6mm	4.50mm
469-475	0.6mm	4.75mm

***Piston shaft length reduced.**

Cat.#	Shaft Dia.	Length
467-375	0.5mm	3.75mm
467-400	0.5mm	4.00mm
467-425	0.5mm	4.25mm
467-450	0.5mm	4.50mm
467-475	0.5mm	4.75mm
467-500	0.5mm	5.00mm
467-600	0.5mm	6.00mm
468-375	0.6mm	3.75mm
468-400	0.6mm	4.00mm
468-425	0.6mm	4.25mm
468-450	0.6mm	4.50mm
468-475	0.6mm	4.75mm
468-500	0.6mm	5.00mm
468-600	0.6mm	6.00mm