


K Kraus Helix

Endoskeletal Ossicular Reconstruction

Innovative K-Helix coil design permits universal adjustability in the OR

Normal anatomy, such as a non-diseased incudo-malleolar joint, is preserved

No Crimping Necessary



Kraus K-Helix piston may be used with a mobile or immobile footplate

Optional cement application enhances prosthesis stability during healing.



Ref #	Description
756-300	K Helix Crown
757-300	K Helix Piston
636-060	.6mm Footplate Shoe

Designed by Eric M. Kraus, M.D., M.S., F.A.C.S. The Ear Center of Greensboro, P.A., Greensboro, NC

8500 Wolf Lake Drive, Suite 110 • Memphis, TN 38133 • USA Toll Free 1.866.472.2363 • www.gracemedical.com

LIT0075 CID2892