











KaWe-polttimot					
	Laite	Tuotenumero	Muu koodi	Jännite	Tyyppi
	Dialight C, Piccolight C, Eurolight C10/ C30/ VET C30, Piccolight F.O. ja F.O.30	99-10525	12.75111.013 (28943)	2.5 V	Vacuum
		99-90105	12.75122.013 (28937)	2.5 V	Krypton
		99-90106	12.75131.013 (28933)	2.5 V	Halogen
		99-90204	12.75231.013 (28922)	3.5 V	Halogen
		99-90228	12.75241.003 (28925)	3.5 V	Xenon
	Eurolight E30/ E50, Piccolight D, Eurolight D30/ E36/ E56	99-13101	12.75112.003 (28944)	2.5 V	Vacuum
		99-13104	12.75132.003 (28946)	2.5 V	Halogen
		99-13103	12.75232.003 (28923)	3.5 V	Halogen
	Eurolight E10	99-13102	12.75113.003 (28947)	2.5 V	Vacuum
	F.O. Laryngoskoopin kädensija	99-90103	12.75121.003 (28935)	2.5 V	Krypton
		99-90127	12.75142.013 (28952)	2.5 V	Xenon
		99-90226	12.75242.013 (28926)	3.5 V	Xenon
		99-90110	12.75123.013 (28941)	2.5 V	Krypton
	Standard-laryngoskoopin kielet 00, 0, 1	99-867	12.75126.003 (28958)	2.5 V	Krypton
	Standard-laryngoskoopin kielet 2, 3, 4, 5	99-85204	12.75127.003 (28959)	2.5 V	Krypton
	F.O. Laryngoskoopin kädensija	99-90106	12.75131.013 (28933)	2.5 V	Halogen
		99-90127	12.75141.003 (28934)	2.5 V	Xenon
		99-90127	12.75142.013 (28952)	2.5 V	Xenon
		99-90204	12.75231.013 (28922)	3.5 V	Halogen
		99-90228	12.75241.003 (28925)	3.5 V	Xenon
		99-90128	12.75143.013 (28954)	2.5 V	Xenon
		99-90227	12.75243.013 (28928)	3.5 V	Xenon
		99-90211	12.75233.013 (28930)	3.5 V	Halogen