



















## Riester-polttimot

	Laite	Tuotenumero	Muu koodi	Jännite	Tyyppi
	Otoskooppi UNI I, speculight	99-114 99-1120 99-1118	10421 10590 11592	2.7 V 2.5 V 3.5 V	Vacuum Halogen Xenon
	Oftalmoskooppi May	99-0045 99-1117 99-1119	10424 10591 10593	2.7 V 2.5 V 3.5 V	Vacuum Halogen Xenon
	Otoskooppi Penscope	99-2060	10472	2.5 V	Vacuum
	Oftalmoskooppi Penscope	99-184	10473	2.7 V	Vacuum
	Otoskoopit: ri-mini, ri-star, ri-scope	99-0600 99-0601 99-1107	10600 10601 10607	2.5 V 3.5 V 3.5 V	Halogen Halogen Xenon
	Operaatio-otoskooppi ri-scope	99-0602 99-0603 99-1109	10602 10603 10609	2.5 V 3.5 V 3.5 V	Halogen Halogen Xenon
	Oftalmoskoopit: ri-mini, ri-star, ri-scope, ri-derma	99-0605 99-0606 99-1108	10605 10606 10608	2.5 V 3.5 V 3.5 V	Halogen Halogen Xenon
		99-1110 99-1112 99-1113	10610 10615 10616	3.5 V 2.5 V 3.5 V	Xenon Halogen Halogen
		99-1111 99-1114 99-1115	10611 10620 10621	3.5 V 2.5 V 3.5 V	Xenon Halogen Halogen
		99-1116	11178	2.5 V	Halogen
	Otsalamppu Clar, Normalux	99-3724	11301	6.0 V	Vacuum
	Otsalamppu ri-clar	99-320	11320	6.0 V	Halogen
	Laryngoskoopin kielet koot 00 – 1	99-169	11378	2.7 V	Vacuum
	Laryngoskoopin kielet koot 2 – 4	99-171	11379	2.7 V	Vacuum
	Laryngoskoopin kielet koot 00 – 1 <i>vanha malli</i>	99-5029	11380	2.7 V	Vacuum
	Laryngoskoopin kielet koot 2 – 4 <i>vanha malli</i>	99-5030	11381	2.7 V	Vacuum

	F.O. laryngoskooppi	99-0630 99-1121 99-1122	11430 11428 11429	2.5 V 2.5 V 3.5 V	Halogen Xenon Xenon
	Laryngoskooppi toplight	99-90106	11431	2.5 V	Halogen