

MEDIVATORS®
ENDO SMARTCAP®
ENDOGATOR®
Lens Cleaning and Irrigation Disposables



MEDIVATORS®
ENDO SMARTCAP®
ENDOGATOR®
Lens Cleaning and Irrigation Disposables

Disposable GI Endoscopy Irrigation Solutions designed with patient safety in mind.

Infection Prevention

- No need for reprocessing
- Reduced risk of cross contamination

Efficiency

- Minimizes reprocessing requirements
- Creates consistent practices

Compatibility

- Use of CO₂ or air insufflation
- Works with all irrigation pumps and sterile water bottles



MEDIVATORS® ENDO SMARTCAP®

Disposable Tubing

- Disposable alternative to the reusable water bottle
- Increases infection prevention measures



MEDIVATORS ENDOGATOR®

Disposable Irrigation Tubing

- Cost effective and completely disposable
- 24 hour use without the need to reprocess
- Reduces infection prevention concerns



MEDIVATORS ENDOGATOR®

Hybrid Tubing

- Combined lens cleaning and irrigation system
- Uses room air or CO₂
- Fits any sterile water bottle



MEDIVATORS ENDOGATOR™

Kit

- One order number
- Innovative lens cleaning/irrigation system
- Option to use room air or CO₂ insufflation

MEDIVATORS® ENDO SMARTCAP® Disposable Tubing

Model Numbers US OUS	Description	Qty/Box
100145 100145U	For OLYMPUS® 140/160/180/190 Series GI Endoscopes (CO ₂ Also Available)	10
100150 100150U	For OLYMPUS 100/130 Series GI Endoscopes (CO ₂ Also Available)	10
100160 100160U	For PENTAX® GI Endoscopes (CO ₂ Also Available)	10
100165 100165U	For FUJIFILM™ GI Endoscopes (CO ₂ Also Available)	10
100141 100141	Adapter for FUJIFILM 530 Series Colonoscopes	50

MEDIVATORS ENDOGATOR® Disposable Irrigation Tubing

Model Number US OUS	Description	Qty/Box
100130 100130U	For ENDO STRATUS®, OLYMPUS OFP, or ENDOSTAT™ II Irrigation Pumps	10
200230 200230U	For OLYMPUS AFU-100, OLYMPUS OFP, ERBE® EIP2 or EGP-100 Irrigation Pumps	10
200230F 200230FU	For EGP-100 Irrigation Pump with FUJIFILM Connector	10
100241 100241	Single Use Auxiliary Water Jet Connector for OLYMPUS 140/160/180/190 series GI Endoscopes	50
100242 100242	Single Use Auxiliary Water Jet Connector for PENTAX GI Endoscopes	50
100130F 100130FU	For ENDO STRATUS, OLYMPUS OFP, or ENDOSTAT™ II Irrigation Pumps for FUJIFILM GI Endoscopes	10
100135 100135	Channel Irrigation Tubing	10

MEDIVATORS ENDOGATOR Hybrid Tubing

Model Number US OUS	Description	Qty/Box
100605 100605S	Hybrid Tubing with CO ₂ connection for OLYMPUS 140/160/180/190 Series GI Endoscopes with ENDO STRATUS, OLYMPUS OFP or ENDOSTAT II Irrigation Pumps with 24 hour auxiliary water jet connector	10
100606 100606S	Hybrid Tubing with CO ₂ connection for OLYMPUS 140/160/180/190 Series GI Endoscopes with OLYMPUS AFU-100, OLYMPUS OFP, EGP-100 or ERBE EIP2 Irrigation Pumps with 24 hour auxiliary water jet connector	10
100608 100608S	Hybrid Tubing with CO ₂ connection for PENTAX GI Endoscopes with EGP-100 or ERBE EIP2 Irrigation Pumps	10
100609 100609S	Hybrid Tubing with CO ₂ connection for OLYMPUS 140/160/180/190 Series GI Endoscopes with ENDO STRATUS, OLYMPUS OFP or ENDOSTAT II Irrigation Pumps	10
100610 100610S	Hybrid Tubing with CO ₂ connection for OLYMPUS 140/160/180/190 Series GI Endoscopes with OLYMPUS AFU-100, OLYMPUS OFP, EGP-100 or ERBE EIP2 Irrigation Pumps	10

MEDIVATORS ENDOGATOR Kit Tubing

Please contact your MEDIVATORS consultant for kit options

ENDO SMARTCAP®, ENDOGATOR® AND ENDO STRATUS® are trademarks and registered trademarks of Medivators Inc.
 PENTAX® is a registered trademark of Hoya Corporation. FUJIFILM™ is a registered trademark of Fujifilm Corporation.
 OLYMPUS® is a registered trademark of Olympus Corporation. ERBE® is a registered trademark of ERBE Elektromedizin GmbH.
 ENDOSTAT™ is a trademark of Boston Scientific Corporation.



Medivators Inc.
 14605 28th Avenue North
 Minneapolis, MN
 55447-4822 USA
 Tel: +1.763.553.3300
 Fax: +1.763.553.3387

Medivators BV
 Sourethweg 11
 6422 PC Heerlen
 The Netherlands
 Tel: +31.45.5.471.471
 Fax: +31.45.5.429.595

Cantel Medical (UK) Limited
 Wolsley House
 Staffordshire Technology Park
 Beaconside, Stafford ST18 0GA
 United Kingdom
 Tel: +44 (0) 1785.782420
 Fax: +44 (0) 1785.782421

Cantel Medical Asia/Pacific Pte. Ltd.
 1A International Business Park
 #05-01 Singapore 609933
 Tel: +65.6227.9698
 Fax: +65.6225.6848

Medivators Inc. Beijing Representative Office
 Room 1801, Floor 18th, Tower A,
 Beijing Marriott Hotel, Office Building
 No. 7, Jianguomen South Avenue,
 Dongcheng District, Beijing 100005 China
 Tel: +86 10 65204039

www.medivators.com

© 2014 Medivators Inc. P/N 50098-246/B