



L-Line SAMURAI

Micro Scissors

KEY FEATURES
OF THE SAMURAI BLADE


- to open the arachnoid mater
- to open the sylvian fissure
- very fine dissection
- extremely sharp
- precise cut

STRAIGHT HANDLE




64.016.41

Blade sideward curved, pointed, blade 12 mm, working length (a) 65 mm, 2 5/8" total length (b) 180 mm, 7"




BAYONET HANDLE




70.017.01

Blade sideward **right curved**, pointed, blade 12 mm, working length (a) 90 mm, 3 1/2" total length (b) 220 mm, 8 3/4"



70.017.11

Blade sideward **left curved**, pointed, blade 12mm, working length (a) 90 mm, 3 1/2" total length (b) 220 mm, 8 3/4"



MEMORY-SHAFT



46.627

Blade sideward curved, pointed, blade 12 mm, working length (a) 90 mm, 3 1/2"




46.637

Blade sideward curved, pointed, blade 12 mm, working length (a) 110 mm, 4 3/8"



- 360° ROTATABLE WORKING ENDS**
- Hand controlled wheel ensures exact positioning of the working end
 - Rotation combined with malleable shaft offers unparalleled instrument positioning



- MEMORY SHAFT (NITINOL)**
- Shapes to coincide with anatomical needs
 - Reverts to original shape after sterilization





HEADQUARTERS

Peter Lazic GmbH

Microsurgical Innovations
Immelmannweg 2
78532 Tuttlingen
GERMANY

Tel. +49 (0) 7461 96643-0
Fax +49 (0) 7461 8745
info@lazic.de
www.lazic.de

Peter Lazic US Inc.

437 Portsmouth Ave
Unit 1B
Greenland, NH 03840
USA

Local: 603 427 8097
Fax: 603 501 0379
support@peterlazic.com
www.peterlazic.com