


SAUGERMED

Surgical Brain Pads

CE 0476

For neurosurgical operations special produced haemostatic featured pads.



Enables safe use

Also, all the pads have an atraumatic featured Trace Rope with smooth touch. Pad locations can easily be traced and identified by these trace ropes. Pls. see figure 1

To avoid any distention Trace rope looped with a special technics*. Pls. see figure 2

All the pads consists of radio-pack blue strap which is X-ray detectable. Pls. see figure 3

To recognize the number of used pads during the operation it has a Special Counter Card 1 to 10. Pls. see figure 4

To prevent the pads from deformation it consists of a protective belt.

Pls. see figure 5

Atraumatic

Produced from hydrofiber structured organic and inorganic materials in nonwoven form Atraumatic and soft.

Unfuzzy and unfeathered, does not leave loofah.


Sterile and dirtiness protected

Gamma sterilized.

Tvyek made of sterile pouch prevents unbounded tear and can only be open from its opening part.


It has an anticontamination bag in the sterile pouch.

There are 16 different dimensions of each pad made of 11 various raw material.




* Patent rights of the special loop protected by the certification of the Turkish Patent Institute as Utility Model.

Surgery & Operation Products




1:1 scale product table (mm.)

Cotton Pads




SAUGERMED Pure Cotton

- Absorbation capacity 6.3 times of its weight.
- Very thin, very light structured.
- Made of 100% cotton raw materials




SAUGERMED Pure Cotton Thick

- Absorbation capacity 8.3 times of its weight
- Smooth , delicate structured.
- Made of 100% cotton raw materials



SAUGERMED Pure Cotton Medium

- Absorbation capacity 7.6 times of its weight
- Smooth , delicate structured.
- Made of 100% cotton raw materials



SAUGERMED Pure Cotton Soft

- Absorbation capacity 11.5 times of its weight
- Very soft, very absorbative, closely structured.
- Made of 100% cotton raw materials.



Viscose Pads




SAUGERMED Thinpad

- Absorbation capacity 7.1 times of its weight.
- Very thin, smooth and silky structured.
- Made of viscose and synthetic additive material.




SAUGERMED Plus

- Absorbation capacity 6 times of its weight .
- Thin, closely structured.
- Made of viscose and synthetic additive material.



SAUGERMED Thincot

- Absorbation capacity 8.5 times of its weight.
- Thin, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




SAUGERMED Laminat

- Absorbation capacity 8.3 times of its weight.
- Standard thickness, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




SAUGERMED Soft

- Absorbation capacity 7.3 times of its weight
- Standard thickness, closely structured, laminated surface.
- Made of viscose and synthetic additive material.



SAUGERMED Ultra Soft

- Absorbation capacity 8.5 times of its weight.
- Smooth , delicate structured.
- Made of viscose and synthetic additive material.



SAUGERMED Cotton Soft

- Absorbation capacity 10 times of its weight.
- Very smooth, delicate structured.
- Made of viscose and synthetic additive material.